


# ProScope HR Software v2.0 - QuickStart

## STEP 1

After downloading, move the ProScope HR v2.x app to your Applications folder.


## STEP 2

On first run, you may need to change your Security & Privacy settings to allow running apps not downloaded from the Mac App Store.




## STEP 3

Change the Save To location to a new file location of your choosing. Optionally, change the Prefix for pic names.



## STEP 4

After plugging in your ProScope HR microscope, choose ProScope HR from the Camera Source menu. Optionally, change the Resolution setting for the camera.



## STEP 5

Capture an image by pressing the capture button on the ProScope HR / HR2 or by clicking the Capture button in the software. Double-click the resulting thumbnail to bring up the pic in the Preview app.

